

6-1-2009

Section: Chemistry

TANTALUM AS AN INDICATOR ELECTRODE IN DIFFERENT POTENTIOMETRIC TITRATIONS

ENAS ATTIA

Chemistry Department, Faculty of Science (Girls), Al- Azhar University, Cairo, Egypt.

Follow this and additional works at: <https://absb.researchcommons.org/journal>



Part of the [Life Sciences Commons](#)

How to Cite This Article

ATTIA, ENAS (2009) "TANTALUM AS AN INDICATOR ELECTRODE IN DIFFERENT POTENTIOMETRIC TITRATIONS," *Al-Azhar Bulletin of Science*: Vol. 20: Iss. 1, Article 7.

DOI: <https://doi.org/10.21608/absb.2009.7545>

This Original Article is brought to you for free and open access by Al-Azhar Bulletin of Science. It has been accepted for inclusion in Al-Azhar Bulletin of Science by an authorized editor of Al-Azhar Bulletin of Science. For more information, please contact kh_Mekheimer@azhar.edu.eg.

TANTALUM AS AN INDICATOR ELECTRODE IN DIFFERENT POTENTIOMETRIC TITRATIONS

ENAS M. ATTIA

Chemistry Department, Faculty of Science (Girls), Al- Azhar University, Cairo, Egypt.

Abstract

Tantalum metal was employed as an indicator electrode in different potentiometric titrations. Silver, platinum and SCE were successfully applied as reference electrodes. A number of classical acid/base reactions were chosen. This involve titration of strong and weak mono- and poly basic acids; HCl, H₂SO₄, H₃PO₄, HCOOH and H₂C₂O₄ with NaOH and NH₄OH. Both forward and backward titrations were carried out. The titration performed until the greatest jump of potential appeared by adding one drop of titrant. Both forward and backward titration curves showed marked potential inflexions which recorded at the theoretically calculated equivalence points. Bimetallic combinations were equally suitable for the determination of equivalence points of acids neutralizing in more than one step, *e.g.*, H₃PO₄ but not with oxalic acid.

The titrations of S₂O₃²⁻ with I₂ solution using the different tantalum systems were equally successful; sharp potential jumps and distinct end points were obtained. Moreover the potential changes at the same direction as that of Pt/SCE system. On the other hand, when ferrous iron reacted with permanganate, a Ta/Pt system gave a potential responsive to Fe⁺² concentrations. The titration curves obtained exhibited regular defined and sharp inflexion points with high potential drop at the equivalence point. This potential drop had a reverse direction compared with that of Pt/SCE system. While on titrating Fe⁺⁺ ions with dichromate solution, the curves had no distinct end points at any of the examined systems. These systems successes also to define the end point in the titration of FeCl₃ against NaOH solution.

Introduction

In a potentiometric titration the endpoint is determined by the use of a pair of electrodes or a combination electrode. The endpoint occurs where there is a maximal rate of change of potential at the endpoint of the titration. In the portion of the curve corresponding to this large change in potential, there is a point at which the curve changes its direction of curvature. This point is an inflection point or break in the curve and ideally it occurs at the equivalence point of the titration [1]. However, in certain cases, there may be a bias in the analysis, and the break may be slightly displaced from the true equivalence point [2]. The electrodes for potentiometric titrations are chosen so that a change in potential of the titration solution, caused by titration of the analyte of interest, is optimally detected. Indicator electrodes qualify

certain requirements depending on the reaction to be tested. In neutralization reactions, the electrode should be both acid and alkali resistant, and respond linearly to variations in pH [3,4]. In redox titrations, the metal oxide layer must possess high electrical conductance. The tantalum electrode and a large number of metal– metal oxide electrodes are commonly used in this type of titrations [5,6]. A potentiometric titration technique has a variety of important applications in our life, [it was used for characterizing vegetable oil quality control \(olive and seed oils\)](#) [7]. Equilibrium studies on the nickel(II) complexes of oxygen, nitrogen and thioether sulfur donor ligands were performed by potentiometric titration [8]. Also, the equilibria of the complexation processes of V^{3+} with aromatic amino acids [9], the protonation and stability constants toward Co^{2+} , Ni^{2+} , Cu^{2+} , Zn^{2+} , Cd^{2+} and Pb^{2+} were determined by potentiometric technique [10]. [Determination of chlorine and dissolved oxygen in waters by iodimetric potentiometric titration using a membrane containing poly vinyl chloride and 2-nitrophenyl octyl ether as plasticizer](#) was determined by means of potentiometric measurements [11]. PVC based liquid membrane electrode coatings is also used in potentiometric detection for a series of substances in cation-exchange systems [12]. This technique has greatest applications in evaluating the detection limit in ion selective electrodes [13-16] and in determining the low concentrations of anionic [17-21] and cationic surfactants [22].

Acid/ base potentiometric titrations have been recently used to determine the dissociation constants of carboxylic groups [23] and to determine the copolymer compositions of *N*-isopropylacrylamide and itaconic acid [24]. The amount of surface active —OH groups and surface charge density on the red mud particles generated from Chinese diaspore bauxite were also evaluated from the acid/basic potentiometric titration data in 0.5 mol/L NaCl solution [25]. The use of Ta/Ta₂O₅ electrode as an indicator electrode in potentiometric acid-base titrations in fused KNO₃ [26] and as a pH indicator electrode in aqueous media [27] was of great interest.

The purpose of the present work is to establish the feasibility of tantalum metal as an indicator electrode with SCE as a reference electrode; and both of Ta/Pt and Ta/Ag bimetallic combinations to be used for potentiometric determination of equivalence points in different types of reactions.

Experiment

Preparation of the indicator-electrode

The tantalum electrode was prepared from a cylindrical rod of spectroscopically pure tantalum metal of 0.096 cm² area supplied by Johnson Matthey, England. A stout copper wire lead was affixed at one end by mechanical jamming. The electrode surface was mechanically polished using 320, 600, 800 and 1000 grid emery papers. Then, it was rubbed against soft cloth to acquire a silvery bright appearance.

Preparation of reference and metal electrodes

The potential was measured relative to saturated calomel electrode (SCE). A silver wire electrode of 0.12cm² area was cleaned to brightness using silver polish on a damp tissue. Platinum sheet with area of 4cm² was cleaned using an aluminum oxide polishing strip. The electrodes were thoroughly rinsed with bi-distilled water.

Preparation for titration

The chemicals required were of analytical reagent grade and were prepared by dilution from a concentrated stock solution and standardized according to recommended procedures [28]. Solutions were stirred on a magnetic stirrer, the electrodes and burette tip were immersed in the solutions. The stirrer speed was adjusted to stir rapidly without splashing or creating a vortex.

The titrant was added in a dosage rate of 0.5 ml/min. The meter readings were taken with each succeeding 0.2 ml increment with making at least from five to ten more additions of titrant before ending the titration. Each system was titrated at least three times to ensure reproducibility. An initial titration was done. This allowed the electrode the opportunity to become equilibrated in the system in which it will be used discarding the results from this initial titration. For the sake of comparison, the titrations were followed simultaneously with Pt/ SCE cell.


In forward titrations, the base was added to 3ml of the acid whereas the reverse was done in the backward titrations. In each experiment, tantalum electrode was left in the solution to be titrated till the establishment of the steady state potential. The oxidation– reduction study involved the titration of I₂ solution with Na₂S₂O₃ and titration of Fe(NH₄)₂(SO₄)₂ against KMnO₄ and K₂Cr₂O₇. With all reagents both

forward (in which oxidant is the titrant) and backward (in which the reducing agent is the titrant) titrations, were carried out. Sodium hydroxide solutions were prepared from carbonate-free NaOH and standardized against oxalic acid using phenolphthalein indicator in order to precipitate ferric iron from FeCl₃ solution. The EMF is plotted against the volume of titrant. A digital potentiometer with a sensitivity of 0.1mV (KEITHLEY, model 175, USA) was used for potential measurements.

Results and Discussion

1. Neutralization reactions

In forward titrations of all acids under test, the potential of tantalum electrode relative to that of SCE, Pt or Ag is somewhat of a negative value at the beginning of the neutralization. However, this potential becomes more negative as the neutralization goes to its completion. The titration curves are characterized by a single potential inflexion at the equivalence point. HCl is a mineral acid which appears to be completely ionized in water according to the reaction:



For potentiometric titration of HCl against NaOH (Figure 1), the stoichiometric point is at $[\text{H}^+] = [\text{OH}] = \frac{1}{2} K_w = 7$ since at this point the solution contains the ions of NaCl which does not change the reaction of pure water [29]. The behavior of H₂SO₄ is similar to that of HCl (Figure 2).




Figure (1): Forward and backward potentiometric titration curves of 0.0799N HCl # 0.0857N NaOH solution.




Figure (2): Forward and backward potentiometric titration curves of 0.0943N H₂SO₄ # 0.0857N NaOH solution.

Forward titration curves of phosphoric acid against NaOH solution using the different electrode systems are characterized by two potential inflexions (Figure 3).






Figure (3): Forward and backward potentiometric titration curves of 0.0914N H₃PO₄ # 0.0857N NaOH solution.

The first inflexion corresponds to one- third of the total replaceable hydrogen. The second appears after two thirds of the replaceable hydrogen has been titrated. The uptake of the titrant along the two steps is equal and corresponds to 1 mole NaOH/ 1 mole H₃PO₄, representing the successive formation of mono- and di- basic salts according to:



The splitting of the two neutralization steps is in accordance with the marked difference between K_1 and K_2 of H₃PO₄, which are amounting to 1.1×10^{-2} and 7.5×10^{-8} respectively [29]. No inflexion

TANTALUM AS AN INDICATOR ELECTRODE IN DIFFERENT S..... 77

is recorded corresponding to the third ionization. This is because the last third of the hydrogen of phosphoric acid is so slightly ionized and the ionization constant is so near that of water ($K_3 = 5.0 \times 10^{-13}$). In back titrations there are no sharp change in potential (as in forward titration) with small additions of acid near the end point.

Despite of oxalic is a dibasic acid with $K_1 = 3.8 \times 10^{-2}$ and $K_2 = 6.1 \times 10^{-5}$ [29] which are neutralized along two well defined steps, only one neutralization step is recorded (Figure 4). This step corresponds to reaction of 2mole NaOH / 1mole acid, which reflects the failure of these combinations to split the two neutralization steps.






Figure (4): Forward and backward potentiometric titration curves of 3ml 0.1N H₂C₂O₄ # 0.0857N NaOH solution.

Formic acid is a weak acid ($K_a = 1.8 \times 10^{-4}$). At the equivalence point, there is a solution of sodium formate of some particular concentration [29]:



However, in case of formic acid the value of potential inflexion is the smallest one corresponding to all other acids (Figure 5). Table 1 illustrates the height of potential jump for all systems used in forward titration descending in the following order:



This is attributed to the considerable difference in the strengths of these acids which follow the same order i.e. $\text{pK} = -7, -2, 1.23, 2.12$ and 3.75 for hydrochloric, sulfuric, oxalic, phosphoric and formic acid respectively [30].




Figure (5): Forward and backward potentiometric titration curves of 0.0714N HCOOH # 0.0857N NaOH solution.

When HCl is titrated with NH_4OH using the different electrode systems, the forward and backward titration curves are characterized by presence of a single potential inflexion at the equivalence point (Figure 6). The potential inflexions at the two cases are smaller than that obtained in the titration of HCl against NaOH solution. This is attributed to the difference in strengths between the two bases.




Figure (6): Forward and backward potentiometric titration curves of 0.0799N HCl # 0.103N NH₄OH solution.

Table 1: comparison of different electrode systems in aqueous acid- base titrations

Type of acid	Type of base	System	Forward				Backward				
			P _i mV	V _{eq} , ml		P _j mV	V _{eq} , ml				
				Ca	Ex		Ca	Ex			
Strong acids	HCl	Ta/SC E	254	2.8	3.2	186	3.2	3.3			
		Ta/Pt	278		3.0				192	3.8	
		Ta/Ag	243		3.2				215	3.3	
	H ₂ SO ₄	NaOH	Ta/SC E	242	3.3	3.5	274	2.7	2.8		
			Ta/Pt	212		3.3				118	3.1
			Ta/Ag	228		3.2				166	3.1
	H ₃ PO ₄	NaOH	Ta/SC E	164, 66	1.4, 3.4	1.5, 3.3	96	2.7	2.6		
			Ta/Pt	137, 34		4				-	-
			Ta/Ag	101, 48		1.3, 3.3				224	2.8
				3							
				1.5, 3.4							

			Ta/SC	184		2.4	113		3.7
	HCl	NH ₄ OH	E	134	2.3	2.4	68	3.9	3.4
			Ta/Pt	174		2.4	130		3.7
			Ta/Ag						
			Ta/SC	172		3.7	126		2.2
	(COOH) ₂	NaOH	E	144	3.5	3.6	76	2.4	2.6
			Ta/Pt	148		3.4	79		2.3
			Ta/Ag						
Weak acids			Ta/SC	182		2.5	113		3.6
	HCOOH	NaOH	E	118	2.5	2.8	183	3.6	4.0
			Ta/Pt	125		2.6	137		3.2
			Ta/Ag						
	P _i , potential jump		V _{eq} , equivalence volume					Ca, calculated	

Ex, experimental

Table (1) illustrates the points of maximum inflexion in the titration curves were nearly identical with the theoretical quantities of the acids or bases necessary for neutralization process.

Analysis of a mixture of acids which differ greatly in their strengths, e.g. formic and hydrochloric acids using Ta/Ag system is illustrated in Figure 7. From the curves, the value of potential jump for titrating 2ml formic acid (solution B) with NaOH solution is 164 mV at equivalence volume of 1.8 ml. The highest potential jump (271mV) is reported when titrating 2ml HCl (solution A) with NaOH at equivalence volume of 2.2ml. This potential jump was decrease to 110mV when adding 4ml formic acid to 2ml HCl (solution C) at equivalence volume of 5.8 ml, and to 150 mV when 2ml formic acid is added to 4ml HCl (solution D) at equivalence volume of 6 ml. The equivalence volumes of solutions C and D is nearly equivalent to the total volume of acidic solution used. Also, the potential jump in solution D is higher than that of solution C due to the increasing amount of stronger acid. Hence in the present study the two acids combined together and produce acidic solution with a new strength depends on the quantity of each acid present in solution. Finally, the system succeeded to detect the end point for a mixture of acids as a uniform solution but failed to analyze the components of the mixture. This behavior is attributed to the fact that, in the aqueous medium, it is impossible to obtain two or more inflections of the titration curve of a mixture of two or more acids or bases until the value of their pK greatly differ from one another [31]. For this analysis to be successful the two acids should differ in strength by at least 10⁵/1 or the titration must be performed in non- aqueous media [31, 32].




Figure (7): Potentiometric titration curves of mixture of 0.0714N HCOOH and 0.0799N HCl acids with 0.0857N NaOH using Ta/Ag system (A: 2ml HCl, B: 2ml formic, C: 2ml HCl + 4ml formic, D: 2ml formic + 4ml HCl).

2- Oxidation reactions

Normal titration curves were obtained when $S_2O_3^{2-}$ is titrated with I_2 solution using the different Ta- systems. The potential starts from negative values and become more positive upon the addition of the oxidizing agent (Figure 8). The potential changes at the equivalence point are close to the standard potential of the concerned redox couples. The behavior obtained for different tantalum systems could be attributed to the fact that the standard potential of I_2/I^- and $S_2O_3^{2-}/S_4O_6^{2-}$ couples, viz. +0.62 and +0.08V coincide with the potential range where oxide formation or H_2 evolution are not feasible. Thus, the results obtained using tantalum electrode may be explained on the same ground. This interpretation agrees with the previous findings of Saleh [33] for molybdenum disilicide electrode. The determining factor in redox reactions is the ratio of the concentration of the oxidized and reduced forms of the ionic species. The potential, E, acquired by the tantalum electrode at 25 °C is given by:

$$E = E^\circ + \frac{0.0591}{n} \log \frac{[Ox]}{[Red]} \quad (5)$$

The potential of the immersed Ta electrode is thus controlled by the ratio of these concentrations. During the oxidation of a reducing agent (forward) or reduction of

an oxidizing agent (backward), the ratio and therefore the potential changes more rapidly in the vicinity of the end point of the reaction.




Figure (8): Forward and backward titration curves of 2.0, 3.5 and 5.0 ml of 0.1N $\text{Na}_2\text{S}_2\text{O}_3$ against I_2 solutions.

Different behavior is noted when strong oxidizing agents like MnO_4^- and $\text{Cr}_2\text{O}_7^{2-}$ were used. Using Pt/SCE cell in the oxidation of Fe^{++} with KMnO_4 and $\text{K}_2\text{Cr}_2\text{O}_7$

TANTALUM AS AN INDICATOR ELECTRODE IN DIFFERENT S..... 83

solutions illustrates a normal titration curves. Concerning forward titrations, tantalum electrode exhibits different potential values in the same direction as Pt/SCE system starting from lower values and becomes more positive as a result of the oxidizing agent addition (except for Ta/Pt system in KMnO_4 which behaves reversibly). However, reverse manner occurs in the backward titrations for the two oxidizing agents using Pt/SCE system and only for KMnO_4 using Ta/SCE system (Figures 9 and 10).

In KMnO_4 solution, the potential of Ta was measured relative to SCE and Ag in forward titration increases rapidly at the beginning reaching its maximum value followed by a steady state potential then another rapid increase in potential occurs. This may be attributed to that in forward titration; the potential of the electrode becomes more positive upon the addition of the oxidizing agent, until it reaches an equal or very close value to that of oxide film formation. Thereafter, the electrode potential is unaffected by further addition of the titrant which is represented by the steady state potential range in the curves. The final increase in potential may be owing to the formation of iron III sulfate hydrate precipitate. The potential continues to increase with increasing the amount of precipitate formed with each addition of KMnO_4 . Testing the possibility for detection of the end point from this plateau, the volume of the end point does not correspond to the highest potential before the steady part of the curves. The middle point of the steady part corresponds to the volumes: 6.25, 9.25, 12.5 ml. for Ta/SCE system and 4.75, 7.5 and 9.5 ml. for Ta/Ag system in titrating 5, 10 and 15 ml. of Fe^{++} solution. These volumes are so far from the calculated volumes which obtained by normal titration methods.

Back titrations using Ta/SCE achieve the right end point. While for Ta/Ag system, the highest potentials (-43, -10 and -29) corresponding to 4, 7 and 11 ml. respectively for titrating 5, 10 and 15 ml of Fe^{++} are not accurate equivalence points.

On the other hand, in Ta/Pt system (forward titrations), the variation of the potential values of the indicator electrode with the quantity of oxidizing agent used follows typically the potentiometric titration curves. At the equivalence point, a sharp drop in the potential value of the indicator electrode occurs. This drop is in the reverse direction compared with that of Pt/SCE system. A reverse behavior is

observed with the back titrations. It was observed that a combination of a valve metal and a noble metal gives good results. This is attributed to the highest extent in electron conductivity possessed by such metal pair [34].




Figure (9): Forward and backward titration curves of 5, 10 and 15 ml. of 0.1N of $\text{Fe}(\text{NH}_4)_2(\text{SO}_4)_2$ solution with 0.1N KMnO_4 solution.

TANTALUM AS AN INDICATOR ELECTRODE IN DIFFERENT S..... 85

Analyzing the forward curves of oxidation of Fe^{++} with $\text{K}_2\text{Cr}_2\text{O}_7$ by testing the volumes corresponded to the highest potential using Ta/SCE, Ta/Pt and Ta/Ag systems demonstrated that they are equal to: 4, 9, 14 ml, 2,8,11 ml and 3, 11, 14 ml for titrating 5, 10, 15 ml of Fe^{++} , respectively (Figure 10). These volumes are far from the right ones which ranging about 5, 10 and 15 ml respectively.




Figure (10): Forward and backward titration curves of 5, 10 and 15 ml. of 0.1N of $\text{Fe}(\text{NH}_4)_2(\text{SO}_4)_2$ solution with 0.1N $\text{K}_2\text{Cr}_2\text{O}_7$ solution.

The main product of tantalum is Ta_2O_5 , which grows adhering to the metal substrate, but ruptures when a critical thickness is reached. This sequence of alternate growth and failure is repeated leading to a laminated scale. The thickening during the short

TANTALUM AS AN INDICATOR ELECTRODE IN DIFFERENT S..... 87

Fe ⁺² /Fe ⁺³	KMnO ₄	5.0	648	-	626	-	742	149	661	-
		10.0	701	-	643	-	654	219	491	-
		15.0	668	-	630	-	609	275	495	-
Fe ⁺² /Fe ⁺³	K ₂ Cr ₂ O ₇	5.0	250	-	-	-	211	-	-	-
		10.0	216	-	-	-	123	-	-	-
		15.0	186	-	-	-	83	-	-	-
S ₂ O ₃ ²⁻ /S ₄ O ₆ ²⁻	I ₂	2.0	142	563	440	566	215	287	93	429
		3.5	151	414	325	497	190	214	104	322
		5.0	141	446	340	494	163	229	82	355

From the above observations, the shape of the titration curves depends upon the type of the oxidizing and reducing agents used as well as on the direction in which the titration technique was carried out. Moreover, the way in which the electrodes respond to [Ox]/ [Red] ratio before and after the equivalence point differ in most cases from that recorded for Pt/calomel electrode.

3- Application of tantalum electrode in potentiometric titration of ferric iron

In this study, tantalum electrode is applied as an indicator electrode in the precipitation of ferric iron from neutral solution of its chloride by the addition of NaOH solutions. The ferric chloride was standardized by reduction with stannous chloride and then titrated with potassium dichromate solution. In case of potentiometric titration of FeCl₃ with NaOH solution, a potential jump occurs at the quantitative precipitation of ferric hydroxide (Figure 11).




Figure (11): Potentiometric titration curves of 5ml 0.0874N FeCl_3 with 0.0857N NaOH.


It was observed that the potential jumps in all systems used are higher than that of the standard Pt/SCE system (Table 3). Also, the experimental equivalence volumes are so near from the calculated one.

Table (3): Comparison of different electrode systems in the precipitation of ferric chloride by the addition of NaOH solution

System	P_j , mV	V_{eq} , ml	
		Ca	Ex
Pt/SCE	168		5.0
Ta/Nb	219		5.0
Ta/SCE	457	5.1	5.0
Ta/Ag	418		5.2
Ta/Pt	231		4.7

P_j , potential jump V_{eq} , equivalence volume Ca, calculated Ex,
experimental

From the Table, the potential jump descending in the following order:



Because of differences in E° potentials, the measurement of potential differences between any two indicator electrodes allowed the determination of equivalence points accompanied by special change in potential jump for each system.

Conclusion

Bimetallic combinations were equally suitable for the determination of equivalence points of mono- basic acids and also with acids neutralizing in more than one step, *e.g.*, H_3PO_4 but not with oxalic acid. These combinations fail to analyze a mixture of formic and hydrochloric acids. In oxidation of $\text{S}_2\text{O}_3^{2-}$ with I_2 solution, the three Ta- systems prove to determine the exact end point. Concerning the oxidation of Fe^{2+} with KMnO_4 , Ta/Pt system can only be used for detecting the end point in forward and backward titrations. While in oxidation of Fe^{2+} with $\text{K}_2\text{Cr}_2\text{O}_7$, none of the three systems can be used for detecting the end point. These systems succeeded also to define the end point in the titration of FeCl_3 against NaOH solution.

Hence a Ta/Ag and Ta/ Pt couple can be used instead of the usual Pt/calomel cell, to indicate the end point. The use of such bimetallic combination enables the elimination of the calomel electrode with its salt bridge so that a more easily to handle assembly results.

References

1. A. HULANICKI, S. G. Ąb, *Encyclopedia of Analytical Science* (2005) 114.
2. E. M. ATTIA, *Al- Azhar Bull. Sci.* 19(2008) 201.
3. L. XIAO-JIE, B. PENG, F. LIU, Y. QIN, *Sensors and Actuators B: Chemical* 125 (2007) 656.
4. A. IGLESIAS, R. LÓPEZ, D. GONDAR, J. ANTELO, S. FIOL, F. ARCE, *Chemosphere* 76 (2009)107.
5. A. BARAKA, A.I. ABDEL-ROHMAN, E.A. EL-TAHER, *Mater. Chem. Phys.* 9 (1983) 447.
6. A.A. ATTIA, *Electro. Chim. Acta* 47(2002) 1241.
7. L. TORTORA, M. STEFANELLI, M. MASTROIANNI, L. LVOVA, C. D. NATALE, A. D'AMICO, D. FILIPPINI, I. LUNDSTRÖM, R. PAOLESSE, *Sensors and Actuators B: Chemical*, In Press, Accepted Manuscript, Available online 28 May 2009.
8. B. KURZAK, A. WOŻNA, A. KAMECKA, K. KURZAK, *Polyhedron* 27 (2008)3053.
9. H. PODSIADŁY, Z. KARWECKA, *Polyhedron* 28(2009)1568.
10. G. AMBROSI, M. FORMICA, V. FUSI, L. GIORGI, E. MACEDI, M. MICHELONI, R. PONTELLINI, *Inorg. Chim. Acta*, 362 (2009) 3709.
11. C. SÁNCHEZ-PEDREÑO, J.A. ORTUÑO, J. HERNÁNDEZ, *Analy. Chim. Acta* 333(1996)107.
12. G. BAZYŁAK, L. J. NAGELS, *J. Chromatography A* 973(2002)85.
13. L. CHEN, H. JU, X. ZENG, X. HE, Z. ZHANG, *Anal. Chim. Acta* 437 (2001) 191.
14. A.M. MAZLOUM, H. DEGHANI, M. JALAYER, H.R. ZARE, *Anal. Sci.* 20 (2004) 1667.
15. M.K. AMINI, A. RAFI, I. MOHAMMADPOOR-BALTORK, *Anal. Lett.* 35 (2002) 1795.
16. X.B. ZHANG, Z.X. HAN, Z.H. FANG, G.L. SHEN, R.Q. YU, *Anal. Chim. Acta* 562 (2006) 210.
17. L. ZELENKA, M. SAK-BOSNAR, N. MAREK, B. KOVACS, *Anal. Lett.* 22 (1989) 2791.
18. M. GERLACHE, Z. SENTŐRURK, J.C. VIRE, J.M. KAUFFMANN, *Anal. Chim. Acta* 349 (1997) 59.

19. S. MARTINEZ-BARRACHINA, J. ALONSO, L. MATIA, R. PRATS, M. D. VALLE, *Anal. Chem.* 71 (1999) 3684.
20. J. SANCHEZ, M. D. VALLE, *Electroanalysis* 13 (2001) 471.
21. R. MATESIC-PUACA, M. SAK-BOSNARB, M. BILICA, B.S. GRABARIC, *Sensors and Actuators B* 106 (2005) 221.
22. R. MATESIC-PUACA, M. SAK-BOSNARB, M. BILICA, B.S. GRABARIC, *Electroanalysis* 16 (2004) 843.
23. K. TØMMERAAS, W. PER-OLOF, *Carbohydrate Polymers* 77(2009)194.
24. C. ERBIL, B. TERLAN, Ö. AKDEMIR, A. T. GÖKÇEÖREN, *J. European Polymer* 45(2009)1728.
25. Z. KUN-YU, AH. HUI-PING, Z. LI-JUAN, C. QI-YUAN, *Transactions of Nonferrous Metals Society of China* 18(2008)1285.
26. A. BARAKA, A.I. ABDEL-RAZIK, A.I. ABDEL-ROHMAN, *Surf. Technol.* 25(1985)31.
27. S. SHAHINE, M.S. EL BASIOUNY, *J. Electroanal. Chem.*108 (1980)271.
28. VOGEL, *Quantitative Inorganic Analysis, The English Language Book Society, Longman* (1978) p. 600.
29. H. H. WILLARD, N. H. FURMAN. *Elementary Quantitative Analysis. Education Manual EM 457.* (1940)117.
30. P. W. ATKINS, *Physical Chemistry Atkins* 6th ed. Oxford (1998) p.934.
31. A. P. KRESHKOV, A. A. YAROSLAVTSEV *Course of Analytical Chemistry vol.2 quantitative analysis,* (1977) p. 296.
32. A. A. SOLIMAN, M. Sc Thesis, Al- Azhar University (1982)
33. R. M. SALEH, Ph.D. Thesis, Cairo University (1970).
34. S. A. DARWISH, R. SALIM, *Microchem.* 1(1973) 670.
35. U. R. EVANS, *The Corrosion and Oxidation of Metals* 2nd ed. Arnold, London (1976) p.34.